

Product datasheet for **GA105113**

Human VDAC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	VDAC3 CRISPRa kit - CRISPR gene activation of human voltage dependent anion channel 3
Symbol:	VDAC3
UniProt ID:	<u>Q9Y277</u>
Synonyms:	HD-VDAC3; VDAC-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)