

Product datasheet for **GA105138**

Human WIPF1 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | WIPF1 CRISPRa kit - CRISPR gene activation of human WAS/WASL interacting protein family member 1 |
| Symbol: | WIPF1 |
| UniProt ID: | O43516 |
| Synonyms: | PRPL-2; WAS2; WASPIP; WIP |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)