

Product datasheet for **GA105139**

Human EIF4H activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF4H CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 4H
Symbol:	EIF4H
UniProt ID:	Q15056
Synonyms:	eIF-4H; WBSCR1; WSCR1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)