

Product datasheet for **GA105140**

Human **CYLN2 (CLIP2)** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CLIP2 CRISPRa kit - CRISPR gene activation of human CAP-Gly domain containing linker protein 2
Symbol:	CLIP2
UniProt ID:	<u>Q9UDT6</u>
Synonyms:	CLIP; CLIP-115; CYLN2; WBSCR3; WBSCR4; WSCR3; WSCR4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)