

## Product datasheet for **GA105149**

### Human WNT3 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | WNT3 CRISPRa kit - CRISPR gene activation of human Wnt family member 3  |
| Symbol:       | WNT3  |
| UniProt ID:   | <u><a href="#">P56703</a></u>   |
| Synonyms:     | INT4; TETAMS  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

#### Product images:



[View online >](#)