

Product datasheet for **GA105186**

Human ZAP70 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZAP70 CRISPRa kit - CRISPR gene activation of human zeta chain of T cell receptor associated protein kinase 70
Symbol:	ZAP70
UniProt ID:	<u>P43403</u>
Synonyms:	SRK; STCD; STD; TZK; ZAP-70
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)