

Product datasheet for **GA105198**

Human ZNF3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF3 CRISPRa kit - CRISPR gene activation of human zinc finger protein 3
Symbol:	ZNF3
UniProt ID:	P17036
Synonyms:	A8-51; HF.12; KOX25; PP838; Zfp113
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)