

Product datasheet for **GA105214**

Human ZNF23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF23 CRISPRa kit – CRISPR gene activation of human zinc finger protein 23
Symbol:	ZNF23
Components:	GA105214G1, ZNF23 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105214G2, ZNF23 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105214G3, ZNF23 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P17027</u>
Synonyms:	KOX16; Zfp612; ZNF359; ZNF612
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

