

## Product datasheet for **GA105232**

### Human ZNF66 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                            |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| <b>Description:</b>  | ZNF66 CRISPRa kit - CRISPR gene activation of human zinc finger protein 66                                                                                                                                                                                                                                 |
| <b>Symbol:</b>       | ZNF66                                                                                                                                                                                                                                                                                                      |
| <b>Components:</b>   | GA105232G1, ZNF66 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA105232G2, ZNF66 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA105232G3, ZNF66 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>Synonyms:</b>     | FLJ16537                                                                                                                                                                                                                                                                                                   |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:

