

Product datasheet for **GA105235**

Human ZNF74 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF74 CRISPRa kit - CRISPR gene activation of human zinc finger protein 74
Symbol:	ZNF74
Components:	GA105235G1, ZNF74 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105235G2, ZNF74 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105235G3, ZNF74 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q16587</u>
Synonyms:	COS52; hZNF7; ZFP520; ZNF520
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)