

## Product datasheet for **GA105237**

### Human ZNF75A activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	ZNF75A CRISPRa kit - CRISPR gene activation of human zinc finger protein 75a
Symbol:	ZNF75A
Components:	GA105237G1, ZNF75A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105237G2, ZNF75A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105237G3, ZNF75A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q96N20</a></u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)