

Product datasheet for **GA105250**

Human MKRN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MKRN3 CRISPRa kit - CRISPR gene activation of human makorin ring finger protein 3
Symbol:	MKRN3
Components:	GA105250G1, MKRN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105250G2, MKRN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105250G3, MKRN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q13064</u>
Synonyms:	CPPB2; D15S9; RNF63; ZFP127; ZNF127
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)