

Product datasheet for **GA105262**

Human Mel18 (PCGF2) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | PCGF2 CRISPRa kit - CRISPR gene activation of human polycomb group ring finger 2 |
| Symbol: | PCGF2 |
| UniProt ID: | P35227 |
| Synonyms: | MEL-18; RNF110; ZNF144 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)