

Product datasheet for **GA105270**

Human ZNF157 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF157 CRISPRa kit - CRISPR gene activation of human zinc finger protein 157
Symbol:	ZNF157
Components:	GA105270G1, ZNF157 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105270G2, ZNF157 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105270G3, ZNF157 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P51786</u>
Synonyms:	HZF22
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

