

Product datasheet for **GA105295**

Human ZNF215 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF215 CRISPRa kit - CRISPR gene activation of human zinc finger protein 215
Symbol:	ZNF215
Components:	GA105295G1, ZNF215 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105295G2, ZNF215 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105295G3, ZNF215 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UL58</u>
Synonyms:	BAZ-2; BAZ2; ZKSCAN11; ZSCAN43
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)