

Product datasheet for **GA105298**

Human ZNF223 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF223 CRISPRa kit – CRISPR gene activation of human zinc finger protein 223
Symbol:	ZNF223
Components:	GA105298G1, ZNF223 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105298G2, ZNF223 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105298G3, ZNF223 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UK11</u>
Synonyms:	Homo sapiens zinc finger protein 223; KRAB A domain; zinc finger protein 223; Zinc finger protein ZNF223
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

