

Product datasheet for **GA105301**

Human ZNF226 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF226 CRISPRa kit - CRISPR gene activation of human zinc finger protein 226
Symbol:	ZNF226
Components:	GA105301G1, ZNF226 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105301G2, ZNF226 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105301G3, ZNF226 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NYT6</u>
Synonyms:	Kruppel-associated box protein; zinc finger protein 226
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)