

Product datasheet for **GA105304**

Human ZNF229 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF229 CRISPRa kit – CRISPR gene activation of human zinc finger protein 229
Symbol:	ZNF229
Components:	GA105304G1, ZNF229 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105304G2, ZNF229 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105304G3, ZNF229 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UJW7</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

