

Product datasheet for **GA105307**

Human ZNF236 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF236 CRISPRa kit - CRISPR gene activation of human zinc finger protein 236
Symbol:	ZNF236
Components:	GA105307G1, ZNF236 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105307G2, ZNF236 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105307G3, ZNF236 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9UL36
Synonyms:	ZNF236A; ZNF236B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

