

Product datasheet for **GA105308**

Human ZnT1 (SLC30A1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC30A1 CRISPRa kit - CRISPR gene activation of human solute carrier family 30 member 1
Symbol:	SLC30A1
Components:	GA105308G1, ZnT1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105308G2, ZnT1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105308G3, ZnT1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9Y6M5
Synonyms:	ZNT1; ZRC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)