

Product datasheet for **GA105337**

Human CXCR4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CXCR4 CRISPRa kit - CRISPR gene activation of human C-X-C motif chemokine receptor 4
Symbol:	CXCR4
UniProt ID:	P61073
Synonyms:	CD184; D2S201E; FB22; HM89; HSY3RR; LAP-3; LAP3; LCR1; LESTR; NPY3R; NPYR; NPYRL; NPYY3R; WHIM
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)