

Product datasheet for **GA105348**

Human Kv beta 1 (KCNAB1) activation kit by CRISPRa

Product data:

| | |
|---------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | KCNAB1 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily A member regulatory beta subunit 1 |
| Symbol: | KCNAB1 |
| UniProt ID: | Q14722 |
| Synonyms: | AKR6A3; hKvb3; hKvBeta3; KCNA1B; KV-BETA-1; Kvb1.3 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)