

## Product datasheet for **GA105391**

### Human SLC25A16 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SLC25A16 CRISPRa kit - CRISPR gene activation of human solute carrier family 25 member 16
Symbol:	SLC25A16
Components:	GA105391G1, SLC25A16 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105391G2, SLC25A16 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105391G3, SLC25A16 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P16260</u>
Synonyms:	D10S105E; GDA; GDC; HGT.1; hML7; ML7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)