

Product datasheet for **GA105404**

Human SSPN activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SSPN CRISPRa kit - CRISPR gene activation of human sarcospan
Symbol:	SSPN
Components:	GA105404G1, SSPN gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105404G2, SSPN gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105404G3, SSPN gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q14714
Synonyms:	DAGA5; KRAG; NSPN; SPN1; SPN2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)