

Product datasheet for **GA105414**

Human CERD4 (DPF3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DPF3 CRISPRa kit - CRISPR gene activation of human double PHD fingers 3
Symbol:	DPF3
UniProt ID:	Q92784
Synonyms:	BAF45C; CERD4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)