

Product datasheet for **GA105421**

Human Gigaxonin (GAN) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GAN CRISPRa kit – CRISPR gene activation of human gigaxonin
Symbol:	GAN
Components:	GA105421G1, Gigaxonin gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105421G2, Gigaxonin gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105421G3, Gigaxonin gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H2C0</u>
Synonyms:	GANI; KLHL16
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

