

Product datasheet for **GA105451**

Human ZRSR2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZRSR2 CRISPRa kit – CRISPR gene activation of human zinc finger CCCH-type, RNA binding motif and serine/arginine rich 2
Symbol:	ZRSR2
Components:	GA105451G1, ZRSR2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105451G2, ZRSR2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105451G3, ZRSR2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q15696
Synonyms:	U2AF1-RS2; U2AF1L2; U2AFIRS2; URP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

