

Product datasheet for **GA105458**

Human HAP40 (F8A1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	F8A1 CRISPRa kit - CRISPR gene activation of human coagulation factor VIII associated 1
Symbol:	F8A1
Components:	GA105458G1, HAP40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105458G2, HAP40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105458G3, HAP40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P23610</u>
Synonyms:	DXS522E; F8A; HAP40
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)