

Product datasheet for **GA105460**

Human TMEM187 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM187 CRISPRa kit – CRISPR gene activation of human transmembrane protein 187
Symbol:	TMEM187
Components:	GA105460G1, TMEM187 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105460G2, TMEM187 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105460G3, TMEM187 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14656</u>
Synonyms:	CXorf12; DXS9878E; ITBA1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

