

Product datasheet for **GA105468**

Human H3FT (HIST3H3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIST3H3 CRISPRa kit - CRISPR gene activation of human histone cluster 3 H3
Symbol:	H3-4
UniProt ID:	Q16695
Synonyms:	H3.4; H3/g; H3FT; H3t
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)