

Product datasheet for **GA105473**

Human TRRAP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRRAP CRISPRa kit - CRISPR gene activation of human transformation/transcription domain associated protein
Symbol:	TRRAP
UniProt ID:	Q9Y4A5
Synonyms:	PAF350/400; PAF400; STAF40; TR-AP; Tra1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)