


## Product datasheet for **GA105496**

### Human HIST1H2AK (HIST1H2AL) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	HIST1H2AL CRISPRa kit - CRISPR gene activation of human histone cluster 1 H2A family member I
Symbol:	H2AC16
Components:	GA105496G1, HIST1H2AK gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105496G2, HIST1H2AK gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105496G3, HIST1H2AK gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<a href="#">P0C0S8</a>
Synonyms:	dj193B12.9; H2A.i; H2A/i; H2AFI
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)