

Product datasheet for **GA105513**

Human H3FA (HIST1H3A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIST1H3A CRISPRa kit - CRISPR gene activation of human histone cluster 1 H3 family member a
Symbol:	H3C1
UniProt ID:	P68431
Synonyms:	H3/A; H3FA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)