

Product datasheet for **GA105515**

Human H3FC (HIST1H3C) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIST1H3C CRISPRa kit - CRISPR gene activation of human histone cluster 1 H3 family member c
Symbol:	H3C3
UniProt ID:	P68431
Synonyms:	H3.1; H3/c; H3FC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)