

Product datasheet for **GA105520**

Human H3FK (HIST1H3H) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIST1H3H CRISPRa kit - CRISPR gene activation of human histone cluster 1 H3 family member h
Symbol:	H3C10
UniProt ID:	P68431
Synonyms:	H3/k; H3F1K; H3FK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)