

Product datasheet for **GA105527**

Human H4C3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | HIST1H4C CRISPRa kit - CRISPR gene activation of human histone cluster 1 H4 family member c |
| Symbol: | H4C3 |
| UniProt ID: | P62805 |
| Synonyms: | dj221C16.1; H4/g; H4FG |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)