

Product datasheet for **GA105546**

Human PIP5K2 beta (PIP4K2B) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PIP4K2B CRISPRa kit - CRISPR gene activation of human phosphatidylinositol-5-phosphate 4-kinase type 2 beta |
| Symbol: | PIP4K2B |
| UniProt ID: | P78356 |
| Synonyms: | PI5P4KB; PIP5K2B; PIP5KIIB; PIP5KIIbeta; PIP5P4KB |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)