

## Product datasheet for **GA105558**

### Human BCAR3 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | BCAR3 CRISPRa kit - CRISPR gene activation of human BCAR3 adaptor protein, NSP family member |
| Symbol:       | BCAR3  |
| UniProt ID:   | <u><a href="#">Q75815</a></u>  |
| Synonyms:     | NSP2; SH2D3B   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                      |

#### Product images:



[View online >](#)