


## Product datasheet for **GA105565**

### Human MST3 (STK24) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | STK24 CRISPRa kit - CRISPR gene activation of human serine/threonine kinase 24 |
| Symbol:       | STK24  |
| UniProt ID:   | <u><a href="#">Q9Y6E0</a></u>  |
| Synonyms:     | HEL-S-95; MST3; MST3B; STE20; STK3   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)        |

#### Product images:



[View online >](#)