


Product datasheet for **GA105576**

Human **DYRK3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	DYRK3 CRISPRa kit - CRISPR gene activation of human dual specificity tyrosine phosphorylation regulated kinase 3
Symbol:	DYRK3
UniProt ID:	O43781
Synonyms:	DYRK5; hYAK3-2; RED; REDK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)