

## Product datasheet for **GA105594**

### Human SPT3 (SUPT3H) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SUPT3H CRISPRa kit - CRISPR gene activation of human SPT3 homolog, SAGA and STAGA complex component
Symbol:	SUPT3H
UniProt ID:	<u><a href="#">Q75486</a></u>
Synonyms:	SPT3; SPT3L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)