

Product datasheet for **GA105619**

Human LAT3 (SLC43A1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC43A1 CRISPRa kit - CRISPR gene activation of human solute carrier family 43 member 1
Symbol:	SLC43A1
UniProt ID:	O75387
Synonyms:	LAT3; PB39; POV1; R00504
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)