

Product datasheet for **GA105621**

Human PI 3 Kinase p55 gamma (PIK3R3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIK3R3 CRISPRa kit - CRISPR gene activation of human phosphoinositide-3-kinase regulatory subunit 3
Symbol:	PIK3R3
UniProt ID:	Q92569
Synonyms:	p55; p55-GAMMA; p55PIK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)