


Product datasheet for **GA105630**

Human Kv beta 2 (KCNAB2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNAB2 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily A regulatory beta subunit 2
Symbol:	KCNAB2
UniProt ID:	Q13303
Synonyms:	AKR6A5; HKvbeta2; HKvbeta2.1; HKvbeta2.2; KCNA2B; KV-BETA-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)