

Product datasheet for **GA105655**

Human PPFIA3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPFIA3 CRISPRa kit - CRISPR gene activation of human PTPRF interacting protein alpha 3
Symbol:	PPFIA3
UniProt ID:	O75145
Synonyms:	LPNA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)