

Product datasheet for **GA105669**

Human Telethonin (TCAP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TCAP CRISPRa kit - CRISPR gene activation of human titin-cap
Symbol:	TCAP
Components:	GA105669G1, Telethonin gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105669G2, Telethonin gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105669G3, Telethonin gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O15273</u>
Synonyms:	CMD1N; CMH25; LGMD2G; T-cap; TELE; telethonin
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)