

## Product datasheet for **GA105675**

### Human kynurenine 3 monooxygenase (KMO) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	KMO CRISPRa kit - CRISPR gene activation of human kynurenine 3-monooxygenase
Symbol:	KMO
UniProt ID:	<a href="#">O15229</a>
Synonyms:	dj317G22.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)