

Product datasheet for **GA105694**

Human FAM193A activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | FAM193A CRISPRa kit - CRISPR gene activation of human family with sequence similarity 193 member A |
| Symbol: | FAM193A |
| Components: | GA105694G1, FAM193A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105694G2, FAM193A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105694G3, FAM193A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>P78312</u> |
| Synonyms: | C4orf8; RES4-22 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)