


Product datasheet for **GA105702**

Human Phosphatidic acid phosphatase type 2B (PLPP3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLPP3 CRISPRa kit - CRISPR gene activation of human phospholipid phosphatase 3
Symbol:	PLPP3
UniProt ID:	O14495
Synonyms:	Dri42; LPP3; PAP2B; PPAP2B; VCIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)