

Product datasheet for **GA105710**

Human PSMG1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMG1 CRISPRa kit - CRISPR gene activation of human proteasome assembly chaperone 1
Symbol:	PSMG1
UniProt ID:	O95456
Synonyms:	C21LRP; DSCR2; LRPC21; PAC-1; PAC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)